

# Protocol for Leek Leaf Homogenization in the Bullet Blender®

The protocol described in this document is for the use of the Bullet Blender® for the homogenization of leek (*Allium ampeloprasum*) leaves. This protocol does not specify a particular buffer - you may choose which is most appropriate for your downstream application (nucleic acid isolation, protein extraction, etc.).

**Materials Required:** leek leaf, Bullet Blender®, homogenization buffer, pipettor, microcentrifuge tubes, and 0.9-2.0mm stainless steel bead blend (part number SSB14B)

## Instructions


1. Cut leaf into long, thin slices of 200mg or less and place each slice into a microcentrifuge tube.
2. Add a volume of beads equal to the mass of tissue. **NOTE:** 100mg  $\cong$  100 $\mu$ L.
3. Add 2 volumes of buffer to the tube for every mass of sample.
4. Close the microcentrifuge tubes and place them into the Bullet Blender®.
5. Set controls for **SPEED 9** and **TIME 4** minutes. Press **Start**.
6. After the run, remove tubes from the instrument.
7. Visually inspect samples. If homogenization is unsatisfactory, run for another two minutes at speed 10.
8. Proceed with your downstream application.

## SAFETY NOTE!!!

**When using a centrifuge to separate your homogenate from the debris and beads, make sure your tubes are balanced.**



Date 05/06/2011 **Before**



**After**



Quasar Instruments, LLC  
4835 Centennial Blvd.  
Colorado Springs, CO 80919