

# Protocol for *D. melanogaster* Adults Homogenization in the Bullet Blender®

The protocol described in this document is for the use of the Bullet Blender® for the homogenization of *Drosophila melanogaster* adults. This protocol does not specify a particular buffer - you may choose which is most appropriate for your downstream application (nucleic acid isolation, protein extraction, etc.).

**Materials Required:** adult *Drosophila*, Bullet Blender®, homogenization buffer, pipettor, microcentrifuge tubes, and [Zirconium Oxide beads \(0.5mm or 1.0mm\)](#) or [Zirconium Silicate beads \(0.5mm\)](#).


## Instructions

1. Place 10-300mg of flies into microcentrifuge tubes.
2. Add a volume of beads equal to the mass of tissue. **NOTE:** 100mg  $\cong$  100 $\mu$ L.
3. Add 0.2mL to 2.0mL buffer (2 volumes of buffer for every volume of sample).
4. Close the microcentrifuge tubes.
5. Place tubes into the Bullet Blender®.
6. Set controls for **SPEED 8** and **TIME 3** minutes.
7. Remove tubes from the instrument.
8. Visually inspect samples, if homogenization is unsatisfactory, run for another two minutes at the **SPEED 10**.
9. Proceed with your downstream application.

## SAFETY NOTE!!!

**When using a centrifuge to separate your homogenate from the debris and beads, make sure your tubes are balanced.**

## TYPICAL RESULTS



**before (flies only)      after (ZrSiO beads)      after (ZrO beads)**

Latest revision: 24 June 2009



Quasar Instruments, LLC  
4835 Centennial Blvd.  
Colorado Springs, CO 80919