

CAPP[®]
D E N M A R K



M a s t e r s o f p r e c i s i o n

CAPP

D E N M A R K

Symphony of Precise Measurements

THE INTRICACY OF THE SCIENCE WORLD IS REMINISCENT OF A MUSIC PERFORMANCE. AS IN MUSIC, A MASTERPIECE PERFORMANCE CALLS FOR ACCURACY IN EACH NOTE PLAYED - IN SCIENCE, THE PRECISION OF MEASUREMENTS USED IN RESEARCH CAN LEAD EITHER TO A GREAT DISCOVERY OR A FRUITLESS QUEST. CAPP RECOGNIZED THIS AT THE EARLIEST STAGES IN ITS DEVELOPMENT AND COMMITTED ITSELF TO ONLY ONE TASK- DESIGNING PIPETTES THAT WILL BE THE MOST PRECISE INSTRUMENTS IN THE INDUSTRY WITHOUT COMPROMISE.

BY INTRODUCING ADJUSTABLE VOLUME PIPETTES THAT CAN ALSO BE USED WITH REMOVABLE FIXED VOLUME KNOBS WE CREATED A PINNACLE OF ACCURACY THAT IS WIDELY USED IN MANY REFERENCE LABORATORIES WORLDWIDE. THERE IS NO LONGER A REASON TO USE MULTIPLE FIXED VOLUME PIPETTES AS A FRAME OF REFERENCE. NOW YOU CAN HAVE ONE CAPP PIPETTE WITH SEPARATE SPECIFIC FIXED VOLUME CONTROLLERS FOR USE IN MOST CRITICAL APPLICATIONS.

THE CONSTRUCTION OF EACH CAPP PIPETTE HAS A SHARED STORY; A TALE OF ARTISTRY OF THE HANDS WHICH GAVE IT SHAPE, OF CAREFUL MATERIAL SELECTION, OF DEDICATED CARE AND PATIENCE. EACH PIPETTE IS INDIVIDUALLY PRODUCED AND TESTED BY A HIGHLY SKILLED MASTER. IT IS THE PASSION OF THESE MASTERS OF PRECISION AT CAPP WHO GIVE EACH PRODUCT THAT VERY SPECIAL CHARACTERISTIC: A SOUL.

AS WE ARE LOOKING INTO THE FUTURE WE ARE FOCUSING ON DELIVERING A FULL LINE OF LIQUID MEASUREMENT INSTRUMENTS OF THE SAME ROBUST DESIGN AND HIGHEST STANDARD OF PRECISION THAT MADE CAPP PIPETTES THE INSTRUMENT OF CHOICE IN MANY CRITICAL LABORATORIES WORLDWIDE. WE INVITE YOU TO EXPERIENCE THE SYMPHONY OF PRECISE MEASUREMENTS AS POWERFUL AS BEETHOVEN'S NINTH SYMPHONY AND AS ENJOYABLE AS MOZART'S MINUETS.

WE CHALLENGE YOU TO BRING CONFIDENCE TO YOUR MEASUREMENTS. BECOME A CONDUCTOR IN CAPP'S ORCHESTRA OF VIRTUOSOS. EXPERIENCE AND EXPRESS SOUND RESULTS IN EVERY DROP.



CAPP

D E N M A R K

CappComfort

SINGLE CHANNEL PIPETTES

comfort.capp.dk

70% THUMB WORKLOAD REDUCTION WITH CAPPAERO™ COMFORT

EXCESSIVE WORKLOAD OF THE THUMB IS A MAJOR CAUSE OF REPETITIVE STRAIN INJURIES AND USER DISCOMFORT IN PIPETTING.

THE UNIQUE DESIGN OF THE CAPPAERO™ COMFORT SENSATIONALLY REDUCES THE THUMB WORKLOAD BY UP TO 70% COMPARED TO OTHER MANUAL PIPETTES. THIS IS ACHIEVED BY THE DESIGN INNOVATION OF THE EXTERNAL TIP EJECTION HANDLE THAT IS OPERATED NOT BY THUMB BUT YOUR WHOLE FIST; LET YOUR OTHER FINGERS TAKE OVER THE LARGE PART OF THE THUMB WORK-LOAD AND DO WHAT FEELS NATURAL WHILE REDUCING THE FORCE NEEDED FOR ASPIRATION AND BLOW-OUT TO A VERY MINIMUM.


WITH THE CAPPAERO™ COMFORT, NOT ONLY IS THE RISK OF HEALTH PROBLEMS RELATED TO REPETITIVE PIPETTING VIRTUALLY ELIMINATED, BUT BY INCREASING THE EFFICIENCY AND PRECISION IN USER WORK, USER FATIGUE IS GREATLY REDUCED WITH ONE ERGONOMIC SOLUTION!









Ordering Info



CappAero Comfort, SoftLine pipettes


CAT. No.	DESCRIPTION	INACCURACY %	IMPRECISION %	COLOR
C10-1A-SL	CAPPAERO PIPETTE, VARIABLE VOL. 0.5-10 UL, COMFORT SOFT	2.50 /1.00	1.80 /0.50	
C20-1A-SL	CAPPAERO PIPETTE, VARIABLE VOL. 2-20 UL, COMFORT SOFT	5.00 /1.00	1.50 /0.30	
C50-1A-SL	CAPPAERO PIPETTE, VARIABLE VOL. 5-50 UL, COMFORT SOFT	2.00 /1.00	1.20 /0.30	
C100-1A-SL	CAPPAERO PIPETTE, VARIABLE VOL. 10-100 UL, COMFORT SOFT	2.00 /0.70	0.80 /0.20	
C200-1A-SL	CAPPAERO PIPETTE, VARIABLE VOL. 20-200 UL, COMFORT SOFT	1.50 /0.60	0.50 /0.15	
C1000-1A-SL	CAPPAERO PIPETTE, VARIABLE VOL. 100-1000 UL, COMFORT SOFT	1.50 /0.80	0.40 /0.15	
C5000-1A-SL	CAPPAERO PIPETTE, VARIABLE VOL. 1-5 ML, COMFORT SOFT	1.20 /0.60	0.30 /0.15	

CappAero Comfort pipettes

CAT. No.	DESCRIPTION	INACCURACY %	IMPRECISION %	COLOR
C10-1A	CAPPAERO PIPETTE, VARIABLE VOL. 0.5-10 UL, COMFORT	2.50 /1.00	1.80 /0.50	
C20-1A	CAPPAERO PIPETTE, VARIABLE VOL. 2-20 UL, COMFORT	5.00 /1.00	1.50 /0.30	
C50-1A	CAPPAERO PIPETTE, VARIABLE VOL. 5-50 UL, COMFORT	2.00 /1.00	1.20 /0.30	
C100-1A	CAPPAERO PIPETTE, VARIABLE VOL. 10-100 UL, COMFORT	2.00 /0.70	0.80 /0.20	
C200-1A	CAPPAERO PIPETTE, VARIABLE VOL. 20-200 UL, COMFORT	1.50 /0.60	0.50 /0.15	
C1000-1A	CAPPAERO PIPETTE, VARIABLE VOL. 100-1000 UL, COMFORT	1.50 /0.80	0.40 /0.15	
C5000-1A	CAPPAERO PIPETTE, VARIABLE VOL. 1-5 ML, COMFORT	1.20 /0.60	0.30 /0.15	
SK-03-XX	THREE PIPETTES FREE OF CHOICE, STAND, RACKS OF TIPS AND TIMER			



Capp ApS
Landbrugsvej 10 · DK 5260
Odense S · Denmark
tlf +45 6613 6140
fax +45 6613 2770



Quasar Instruments, LLC | 4735 Centennial Blvd, Colorado Springs, CO 80919
Phone: (719) 535-8000; Toll Free (877) 535-8001; Fax: (719) 547-2404
www.QuasarInstruments.com | sales@qi-us.com