

CAPP<sup>®</sup>  
D E N M A R K



*M a s t e r s   o f   p r e c i s i o n*

# CAPP

D E N M A R K

## *Symphony of Precise Measurements*

THE INTRICACY OF THE SCIENCE WORLD IS REMINISCENT OF A MUSIC PERFORMANCE. AS IN MUSIC, A MASTERPIECE PERFORMANCE CALLS FOR ACCURACY IN EACH NOTE PLAYED - IN SCIENCE, THE PRECISION OF MEASUREMENTS USED IN RESEARCH CAN LEAD EITHER TO A GREAT DISCOVERY OR A FRUITLESS QUEST. CAPP RECOGNIZED THIS AT THE EARLIEST STAGES IN ITS DEVELOPMENT AND COMMITTED ITSELF TO ONLY ONE TASK- DESIGNING PIPETTES THAT WILL BE THE MOST PRECISE INSTRUMENTS IN THE INDUSTRY WITHOUT COMPROMISE.

BY INTRODUCING ADJUSTABLE VOLUME PIPETTES THAT CAN ALSO BE USED WITH REMOVABLE FIXED VOLUME KNOBS WE CREATED A PINNACLE OF ACCURACY THAT IS WIDELY USED IN MANY REFERENCE LABORATORIES WORLDWIDE. THERE IS NO LONGER A REASON TO USE MULTIPLE FIXED VOLUME PIPETTES AS A FRAME OF REFERENCE. NOW YOU CAN HAVE ONE CAPP PIPETTE WITH SEPARATE SPECIFIC FIXED VOLUME CONTROLLERS FOR USE IN MOST CRITICAL APPLICATIONS.

THE CONSTRUCTION OF EACH CAPP PIPETTE HAS A SHARED STORY; A TALE OF ARTISTRY OF THE HANDS WHICH GAVE IT SHAPE, OF CAREFUL MATERIAL SELECTION, OF DEDICATED CARE AND PATIENCE. EACH PIPETTE IS INDIVIDUALLY PRODUCED AND TESTED BY A HIGHLY SKILLED MASTER. IT IS THE PASSION OF THESE MASTERS OF PRECISION AT CAPP WHO GIVE EACH PRODUCT THAT VERY SPECIAL CHARACTERISTIC: A SOUL.

AS WE ARE LOOKING INTO THE FUTURE WE ARE FOCUSING ON DELIVERING A FULL LINE OF LIQUID MEASUREMENT INSTRUMENTS OF THE SAME ROBUST DESIGN AND HIGHEST STANDARD OF PRECISION THAT MADE CAPP PIPETTES THE INSTRUMENT OF CHOICE IN MANY CRITICAL LABORATORIES WORLDWIDE. WE INVITE YOU TO EXPERIENCE THE SYMPHONY OF PRECISE MEASUREMENTS AS POWERFUL AS BEETHOVEN'S NINTH SYMPHONY AND AS ENJOYABLE AS MOZART'S MINUETS.

WE CHALLENGE YOU TO BRING CONFIDENCE TO YOUR MEASUREMENTS. BECOME A CONDUCTOR IN CAPP'S ORCHESTRA OF VIRTUOSOS. EXPERIENCE AND EXPRESS SOUND RESULTS IN EVERY DROP.



# CappFixed

FIXED VOLUME PIPETTES

fixed.capp.dk

THERE IS NO LONGER A REASON TO USE MULTIPLE FIXED VOLUME PIPETTES AS A FRAME OF REFERENCE. NOW YOU CAN HAVE ONE CAPP PIPETTE WITH SEPARATE SPECIFIC FIXED VOLUME CONTROLLERS FOR USE IN YOUR MOST CRITICAL APPLICATIONS. THIS UNIQUE DESIGN SUBSTANTIALLY REDUCES INVESTMENT COST AND BENCH SPACE. IT ALSO ALLOWS TO HAVE A SET OF SEPARATE VOLUME CONTROLLERS PRE-CALIBRATED TO ACCURATELY PERFORM WITH DIFFERENT VISCOSITY LIQUIDS. EACH FIXED VOLUME PIPETTE CAN BE CONVERTED INTO VARIABLE VOLUME PIPETTE BY ADDING AN ADJUSTABLE VOLUME CONTROLLER.



## Ordering Info

### Capp ecopipette

CAT. No.	DESCRIPTION	INACCURACY %	IMPRECISION %	COLOR
C02-1S-001	CAPP PIPETTE, FIXED VOLUME 1 UL , ECOPIPETTE	3.10	1.60	
C10-1-005	CAPP PIPETTE, FIXED VOLUME 5 UL, ECOPIPETTE	1.20	0.64	
C10-1-010	CAPP PIPETTE, FIXED VOLUME 10 UL, ECOPIPETTE	1.00	0.50	
C50-1-020	CAPP PIPETTE, FIXED VOLUME 20 UL, ECOPIPETTE	1.00	0.45	
C50-1-025	CAPP PIPETTE, FIXED VOLUME 25 UL, ECOPIPETTE	0.80	0.35	
C100-1-030	CAPP PIPETTE, FIXED VOLUME 30 UL, ECOPIPETTE	1.16	0.40	
C100-1-050	CAPP PIPETTE, FIXED VOLUME 50 UL, ECOPIPETTE	0.90	0.30	
C200-1-100	CAPP PIPETTE, FIXED VOLUME 100 UL , ECOPIPETTE	0.80	0.20	
C200-1-200	CAPP PIPETTE, FIXED VOLUME 200 UL, ECOPIPETTE	0.60	0.15	
C1000-1-500	CAPP PIPETTE, FIXED VOLUME 500 UL, ECOPIPETTE	0.80	0.15	
C1000-1-1000	CAPP PIPETTE, FIXED VOLUME 1000 UL, ECOPIPETTE	0.80	0.15	
C5000-1A-5000	CAPP PIPETTE, FIXED VOLUME 5000 UL, COMFORT	0.60	0.15	
Cxxx-1-yyy	"OTHER FIXED VOL. PIPETTES(NOTE: xxx=MAX PIPETTE VOLUME; yyy=DESIRED VOLUME)"			



Capp ApS  
Landbrugsvej 10 · DK 5260  
Odense S · Denmark  
tlf +45 6613 6140  
fax +45 6613 2770



Quasar Instruments, LLC | 4735 Centennial Blvd, Colorado Springs, CO 80919  
Phone: (719) 535-8000; Toll Free (877) 535-8001; Fax: (719) 547-2404  
[www.QuasarInstruments.com](http://www.QuasarInstruments.com) | [sales@qi-us.com](mailto:sales@qi-us.com)