

CAPP<sup>®</sup>  
D E N M A R K



*M a s t e r s   o f   p r e c i s i o n*

# CAPP

D E N M A R K

## *Symphony of Precise Measurements*

THE INTRICACY OF THE SCIENCE WORLD IS REMINISCENT OF A MUSIC PERFORMANCE. AS IN MUSIC, A MASTERPIECE PERFORMANCE CALLS FOR ACCURACY IN EACH NOTE PLAYED - IN SCIENCE, THE PRECISION OF MEASUREMENTS USED IN RESEARCH CAN LEAD EITHER TO A GREAT DISCOVERY OR A FRUITLESS QUEST. CAPP RECOGNIZED THIS AT THE EARLIEST STAGES IN ITS DEVELOPMENT AND COMMITTED ITSELF TO ONLY ONE TASK- DESIGNING PIPETTES THAT WILL BE THE MOST PRECISE INSTRUMENTS IN THE INDUSTRY WITHOUT COMPROMISE.

BY INTRODUCING ADJUSTABLE VOLUME PIPETTES THAT CAN ALSO BE USED WITH REMOVABLE FIXED VOLUME KNOBS WE CREATED A PINNACLE OF ACCURACY THAT IS WIDELY USED IN MANY REFERENCE LABORATORIES WORLDWIDE. THERE IS NO LONGER A REASON TO USE MULTIPLE FIXED VOLUME PIPETTES AS A FRAME OF REFERENCE. NOW YOU CAN HAVE ONE CAPP PIPETTE WITH SEPARATE SPECIFIC FIXED VOLUME CONTROLLERS FOR USE IN MOST CRITICAL APPLICATIONS.

THE CONSTRUCTION OF EACH CAPP PIPETTE HAS A SHARED STORY; A TALE OF ARTISTRY OF THE HANDS WHICH GAVE IT SHAPE, OF CAREFUL MATERIAL SELECTION, OF DEDICATED CARE AND PATIENCE. EACH PIPETTE IS INDIVIDUALLY PRODUCED AND TESTED BY A HIGHLY SKILLED MASTER. IT IS THE PASSION OF THESE MASTERS OF PRECISION AT CAPP WHO GIVE EACH PRODUCT THAT VERY SPECIAL CHARACTERISTIC: A SOUL.

AS WE ARE LOOKING INTO THE FUTURE WE ARE FOCUSING ON DELIVERING A FULL LINE OF LIQUID MEASUREMENT INSTRUMENTS OF THE SAME ROBUST DESIGN AND HIGHEST STANDARD OF PRECISION THAT MADE CAPP PIPETTES THE INSTRUMENT OF CHOICE IN MANY CRITICAL LABORATORIES WORLDWIDE. WE INVITE YOU TO EXPERIENCE THE SYMPHONY OF PRECISE MEASUREMENTS AS POWERFUL AS BEETHOVEN'S NINTH SYMPHONY AND AS ENJOYABLE AS MOZART'S MINUETS.

WE CHALLENGE YOU TO BRING CONFIDENCE TO YOUR MEASUREMENTS. BECOME A CONDUCTOR IN CAPP'S ORCHESTRA OF VIRTUOSOS. EXPERIENCE AND EXPRESS SOUND RESULTS IN EVERY DROP.



# CAPP

D E N M A R K

## Microbiology

SINGLECHANNEL PIPETTES

[microbiology.capp.dk](http://microbiology.capp.dk)

**CAPP MICROBIOLOGY PIPETTES:** IDEAL FOR MICROBIOLOGY TESTS IN THE FOOD INDUSTRY, THIS PIPETTE IS SUPPLIED WITH THREE FIXED VOLUME KNOBS AT 100UL, 900UL AND 1000UL AND HAS A MODIFIED TIP CONE SUITED FOR USE WITH A STERILE EXTENDED LENGTH STRAW, ALLOWING FURTHER REACH AND PIPETTING OF VISCOUS LIQUIDS SUCH AS MILK, CREAM, YOGHURT AND COSMETICS. IT IS AN IDEAL INSTRUMENTS FOR PROTOCOLS INVOLVING DILUTION OF 1:10.

### MAIN FEATURES:

- EXTENDED REACH•
- ROTATION FREE VOLUME SELECTION•
- PRE-SET CALIBRATED VOLUME CONTROLLERS•
- FULLY AUTOCLAVABLE•

## Ordering Info

### CappAero Microbiology pipettes (for use with sterile straws)

CAT. No.	DESCRIPTION	INACCURACY %	IMPRECISION %
C 1000-1BS	CAPPAERO MICROBIOLOGY PIPETTE, WITH 100UL, 900UL AND 1000UL VOLUME KNOBS	3.00 /1.60	0.80 /0.30
4170807	PRE-STERILE STRAWS FOR CAPPAERO MICROBIOLOGY PIPETTE, 40X25 PCS.		



Capp ApS  
Landbrugsvej 10 · DK 5260  
Odense S · Denmark  
tlf +45 6613 6140  
fax +45 6613 2770



Quasar Instruments, LLC | 4735 Centennial Blvd, Colorado Springs, CO 80919  
Phone: (719) 535-8000; Toll Free (877) 535-8001; Fax: (719) 547-2404  
[www.QuasarInstruments.com](http://www.QuasarInstruments.com) | [sales@qi-us.com](mailto:sales@qi-us.com)