

CAPP[®]
D E N M A R K



M a s t e r s o f p r e c i s i o n

CAPP

D E N M A R K

Symphony of Precise Measurements

THE INTRICACY OF THE SCIENCE WORLD IS REMINISCENT OF A MUSIC PERFORMANCE. AS IN MUSIC, A MASTERPIECE PERFORMANCE CALLS FOR ACCURACY IN EACH NOTE PLAYED - IN SCIENCE, THE PRECISION OF MEASUREMENTS USED IN RESEARCH CAN LEAD EITHER TO A GREAT DISCOVERY OR A FRUITLESS QUEST. CAPP RECOGNIZED THIS AT THE EARLIEST STAGES IN ITS DEVELOPMENT AND COMMITTED ITSELF TO ONLY ONE TASK- DESIGNING PIPETTES THAT WILL BE THE MOST PRECISE INSTRUMENTS IN THE INDUSTRY WITHOUT COMPROMISE.

BY INTRODUCING ADJUSTABLE VOLUME PIPETTES THAT CAN ALSO BE USED WITH REMOVABLE FIXED VOLUME KNOBS WE CREATED A PINNACLE OF ACCURACY THAT IS WIDELY USED IN MANY REFERENCE LABORATORIES WORLDWIDE. THERE IS NO LONGER A REASON TO USE MULTIPLE FIXED VOLUME PIPETTES AS A FRAME OF REFERENCE. NOW YOU CAN HAVE ONE CAPP PIPETTE WITH SEPARATE SPECIFIC FIXED VOLUME CONTROLLERS FOR USE IN MOST CRITICAL APPLICATIONS.

THE CONSTRUCTION OF EACH CAPP PIPETTE HAS A SHARED STORY; A TALE OF ARTISTRY OF THE HANDS WHICH GAVE IT SHAPE, OF CAREFUL MATERIAL SELECTION, OF DEDICATED CARE AND PATIENCE. EACH PIPETTE IS INDIVIDUALLY PRODUCED AND TESTED BY A HIGHLY SKILLED MASTER. IT IS THE PASSION OF THESE MASTERS OF PRECISION AT CAPP WHO GIVE EACH PRODUCT THAT VERY SPECIAL CHARACTERISTIC: A SOUL.

AS WE ARE LOOKING INTO THE FUTURE WE ARE FOCUSING ON DELIVERING A FULL LINE OF LIQUID MEASUREMENT INSTRUMENTS OF THE SAME ROBUST DESIGN AND HIGHEST STANDARD OF PRECISION THAT MADE CAPP PIPETTES THE INSTRUMENT OF CHOICE IN MANY CRITICAL LABORATORIES WORLDWIDE. WE INVITE YOU TO EXPERIENCE THE SYMPHONY OF PRECISE MEASUREMENTS AS POWERFUL AS BEETHOVEN'S NINTH SYMPHONY AND AS ENJOYABLE AS MOZART'S MINUETS.

WE CHALLENGE YOU TO BRING CONFIDENCE TO YOUR MEASUREMENTS. BECOME A CONDUCTOR IN CAPP'S ORCHESTRA OF VIRTUOSOS. EXPERIENCE AND EXPRESS SOUND RESULTS IN EVERY DROP.



 **ecopipette**

SINGLECHANNEL PIPETTE

ecopipette.capp.dk

ecopipette Innovation

- THIS PIPETTE IS BUILT TO LAST. MAKING NO COMPROMISE ON QUALITY IN PRODUCT DEVELOPMENT, RAW MATERIALS USED AND MANUFACTURING PROCESSES WE DEVELOPED THE GREENEST PIPETTE AVAILABLE ON THE MARKET. THE **ECOPIPETTE** IS NOT ONLY PACKAGED USING RECYCLABLE, BIODEGRADABLE MATERIALS, BUT CONSTRUCTED FROM SOME OF THE MOST RENEWABLE DURABLE RESOURCES AVAILABLE TO ENSURE AN ADVANCED LIFE SPAN AND GREENEST AFTER SERVICE RECYCLING.
- CAPP ECOTRADE™ PROGRAM: RETURN ANY BRAND OF OLD PIPETTES IN EXCHANGE FOR DISCOUNTED NEW ECOPIPETTES™. AT CAPP APS, WE WILL ENSURE AN ENVIRONMENTALLY SAFE WASTE DISCHARGE OF ALL OLD PIPETTES. ASK YOUR FAVOURITE AUTHORIZED CAPP DISTRIBUTOR FOR MORE INFORMATION ABOUT OUR ECOTRADE™ PROGRAMME.

CAPP **ECOPIPETTE**™ UNIQUE ADVANTAGES:

- LOCKABLE VOLUME CONTROLLER REQUIRING 50% FEWER TURNS
- DISPLAY IS FULLY VISIBLE DURING HAND OPERATIONS
- REDUCED FORCE REQUIREMENT FOR ASPIRATION AND BLOW-OUT
- THE ONLY PIPETTE WITH FIXED AND VARIABLE EXCHANGEABLE COUNTER MECHANISM
- LARGE PAD FOR EASE OF TIP EJECTION
- FULLY AUTOCLAVABLE PIPETTE BODY (NO DISASSEMBLY REQUIRED)
- DOUBLE O-RING FOR OPTIMAL TIP SEAL
- COMPATIBLE WITH ALL MAJOR PIPETTE TIP BRANDS
- ROBUST DESIGN INCLUDES METAL TIP BASE
- LOWEST INCIDENCE OF BREAKAGE RESULTING FROM ACCIDENTAL DROPPING OR TIPS MOUNTING.






Ordering Info

CAT. No.	DESCRIPTION	INACCURACY %	IMPRECISION %	COLOR
C02-1S	CAPP PIPETTE, VARIABLE VOL. 0.2-2 UL, ECOPIPETTE	12.00 /2.00	6.00 /1.00	
C10-1	CAPP PIPETTE, VARIABLE VOL. 0.5-10 UL, ECOPIPETTE	2.50 /1.00	1.80 /0.50	
C20-1	CAPP PIPETTE, VARIABLE VOL. 2-20 UL, ECOPIPETTE	5.00 /1.00	1.50 /0.30	
C50-1	CAPP PIPETTE, VARIABLE VOL. 5-50 UL, ECOPIPETTE	2.00 /1.00	1.20 /0.30	
C100-1	CAPP PIPETTE, VARIABLE VOL. 10-100 UL, ECOPIPETTE	2.00 /0.70	0.80 /0.20	
C200-1	CAPP PIPETTE, VARIABLE VOL. 20-200 UL, ECOPIPETTE	1.50 /0.60	0.50 /0.15	
C1000-1	CAPP PIPETTE, VARIABLE VOL. 100-1000 UL, ECOPIPETTE	1.50 /0.80	0.40 /0.15	
C5000-1	CAPP PIPETTE, VARIABLE VOL. 1000-5000 UL, ECOPIPETTE	1.20 /0.60	0.30 /0.15	
SKE-03-XX	THREE PIPETTES FREE OF CHOICE, STAND, RACKS OF TIPS AND TIMER			



Capp ApS
Landbrugsvej 10 · DK 5260
Odense S · Denmark
tlf +45 6613 6140
fax +45 6613 2770



Quasar Instruments, LLC | 4735 Centennial Blvd, Colorado Springs, CO 80919
Phone: (719) 535-8000; Toll Free (877) 535-8001; Fax: (719) 547-2404
www.QuasarInstruments.com | sales@qi-us.com