



**Your Source Of Laboratory  
Plasticware.**

# Mini Centrifuge



Mini Centrifuges are ideal for quick spin downs and microfiltration

Each unit has an on/off switch allowing the user to stop the motion of the rotor before opening the centrifuge lid.

- Supplied with rotors for 1.5/ 2.0ml tubes and 0.2ml PCR strips.
- Noisy  $\leq 48$ db.
- 6000 rpm / 2000 x g.
- On / Off switch.
- Durable lid hinge.

## order information

Product No.	Description	Color
E600328	Mini Centrifuges 110V	Green
E610328-B	Mini Centrifuges 220V	Blue

# Extragen

[www.extragen-web.com](http://www.extragen-web.com)



Quasar Instruments, LLC | 4735 Centennial Blvd, Colorado Springs, CO 80919  
Phone: (719) 535-8000; Toll Free (877) 535-8001; Fax (719) 574-2404  
[www.QuasarInstruments.com](http://www.QuasarInstruments.com) | [sales@qi-us.com](mailto:sales@qi-us.com)