

Fantastic perspective for small budgets


MSL4000 - stereo microscopes

MSL4000-series stereo microscopes offer an optimal value for money. Thanks to the wide range of accessories and different eyepieces, these microscopes are suitable for a wide range of applications.

All microscopes have a 45° inclined viewer piece, an eye-distance adjustment and dioptre adjustment. The metal casing is sturdy and durable. To allow you the option for working anywhere without having to depend on any external power supply, some of the MSL microscopes have a battery providing a user-friendly 25 hours of power.

Power supply: 100 - 240 V, 50 - 60 Hz.

A wide range of accessories is available for all models.



MSL4000-10/30-S



MSL4000-10/30-IL-TL

IL Incident light
 TL Transmitted light
 10/30 10x/30x magnification
 20/40 20x/40x magnification
 S Swivelling arm

	Optical equipment	Illumination	Special features
MSL4000-10/30-IL-TL	10x wide-field eyepieces 1x and 3x objectives 10x and 30x magnification	LED incident and transmitted light	Battery life: 25 hours
MSL4000-10/30-IL-S	10x wide-field eyepieces 1x and 3x objectives 10x and 30x magnification	LED incident lighting	Swivelling arm
MSL4000-10/30-S	10x wide-field eyepieces 1x and 3x objectives 10x and 30x magnification		Swivelling arm
MSL4000-20/40-IL-TL	10x wide-field eyepieces 2x and 4x objectives 20x and 40x magnification	LED incident and transmitted light	Battery life: 25 hours
MSL4000-20/40-IL-S	10x wide-field eyepieces 2x and 4x objectives 20x and 40x magnification	LED incident lighting	Swivelling arm
MSL4000-20/40-S	10x wide-field eyepieces 2x and 4x objectives 20x and 40x magnification		Swivelling arm